

Non-Resorbable Aesthetic Injectable Filler for Correction of Nasolabial Folds

In U.S. clinical trials, Artefill's unique composition of polymethylmethacrylate microspheres (20% of implant volume) and purified bovine collagen (80% of implant volume, containing 0.3% lidocaine) provided safe and effective soft tissue augmentation for the correction of nasolabial folds.¹

arte fill®
THE LASTING SOLUTION

Key Benefits

- > A new category of microsphere-enhanced collagen filler for the correction of nasolabial folds
- > Excellent safety profile¹
- > PMMA microspheres have been optimized over 3 generations
- > Unlike temporary dermal fillers, Artefill microspheres create a support matrix for enduring wrinkle correction

Ease of Treatment

- > Contains 0.3% lidocaine, reducing need for nerve blocks and increasing patient comfort
- > Supplied in pre-filled syringes (no reconstitution necessary) that are ready for immediate use
- > Nonporous, proprietary syringes with 0.8 cc fills with optimized flow characteristics
- > Artefill is injected into the dermal/subdermal junction using the tunneling technique. This injection technique is similar to the technique used with many temporary dermal fillers

Proven Results

- > FDA approval based on controlled, randomized, double-masked, prospective, multi-center U.S. clinical trial
- > Patients treated with Artefill (n=251) showed sustained aesthetic improvement throughout the primary study endpoint¹. The collagen control group exhibited no lasting improvements after 6 months
- > Safety profiles for Artefill and the collagen control group were comparable. Artefill patients were evaluated one year after treatment and demonstrated minimal side effects¹

Before Artefill




Actual Artefill Patient

After 6 Months*



* Primary endpoint

After 1 Year**



**No 12-month control per protocol

Specifications

Artefill Treatment Kits:

Artefill Treatment Kit:
AF0508
5 syringes x 0.8 cc
Each treatment kit is packaged with ten 26g x 5/8 inch needles.

Artefill Skin Test Kit:
ST0501
5 syringes
Skin test kit packaged with five 30g x 1/2 inch needles.

Artefill contains:

20% Microspheres
Polymethylmethacrylate
(PMMA)

80% Purified Collagen Gel
Consisting of:

Bovine Collagen	3.5%
Phosphate Buffer	2.7%
Sodium Chloride	0.9%
Lidocaine Hydrochloride	0.3%
Water for Injection	92.6%

Notices:

Sterile, Non-Pyrogenic

Store refrigerated at
35°F to 46°F (2°C to 8°C)

Do not freeze

Single use only

Federal Law restricts this
device to physician use only.

Note: Artefill Skin Test required. Please refer to complete Artefill Instructions for Use for detailed discussion of safety and efficacy.

Ordering Information

To place an order contact Suneva Medical
In the United States call toll free 888.ARTEFILL (888.278.3345)
Orders may also be sent by fax to 858.550.9997

LEGALLY OWNED AND DISTRIBUTED BY



5870 Pacific Center Blvd.
San Diego, CA 92121
United States of America

888.ARTEFILL (278.3345)

858.550.9999 | phone
858.550.9997 | fax

sunevamedical.com

Reference

1. Artefill Instructions For Use. San Diego, CA: Suneva Medical, Inc.; 2009